

SF100C Single Bed Room Station IL362 Specification Sheet

Tek-CARE NC110 Nurse Call System • Rev. 16 - 4/2022

Features

- Simple, push-button calling using plug-in call cord
- Call-placed indicator and cancel button
 - Call automatically registers if call cord is accidentally removed
- Corridor zone light capability

Description

The SF100C Single Bed Room Station includes a jack allowing the attachment of an SF301-series Call Cord. A call is registered at the nurse station, corridor light and corridor zone light (if installed) when the call cord button is used, or if the call cord is accidentally removed from the jack on the panel. The station includes a call-placed indicator that lights when a call is placed and a cancel button mounted on an attractive, flame-retardant (UL[®] 94V-0) plastic faceplate.



SF100C Single Bed Room Station

Required Components

Tek-CARE NC110 Nurse Call System
SF301-series Call Cord

Listings



UL 1069 Listed; UL[®] File Number E71599.

Specifications

Dimensions:	4.5"H × 2.75"W × 1.75"D (114 mm × 69 mm × 44 mm)
Connections:	Plug-in connectors
Construction:	Flame retardant (UL [®] 94V-0) plastic faceplate. Switches and electronic components are circuit board mounted
Housing:	Single-gang ring, or single-gang ring with dual-gang box