

SF301G Call Cord IL434 Specification Sheet

Tek-CARE Nurse Call Systems • Rev. 22- 2/2025

Features

- " Sensitive touch activation
- " Oxygen tent safe
- " Ideal for geriatric use
- Momentary, air-activated switch
- Adjustable bed clip(s)
- UL[®] 1069 Listed

Description

The TekTone[®] SF301G Call Cord is for use with Tek-CARE[®] nurse call systems and wireless emergency call systems where 1/4" phono jacks are used.

The SF301G Call Cord is specifically designed for use in geriatric and oxygen environments, but may also be used in general applications. To place a call, the patient applies slight pressure to the squeeze bulb using any part of the body. The force of air through the vinyl tubing closes a switch in the plug, thus placing a call. An adjustable bed clip is supplied for the SF301G and is attached to the call cord.

Required Components

Tek-CARE[®]3000 Nurse Call System
with SF121 or SF122 Patient Stations
or
Tek-CARE[®]400 GEN3 Nurse Call System
with SF121 or SF122 Patient Stations
or
Tek-CARE[®]500 Wireless Nurse Call System
with SF521UL or SF523UL Room Stations
or
Tek-CARE[®]120 Nurse Call System
with SF121 or SF122 Room Stations
or
Tek-CARE[®]300II Nurse Call System
with SF380A Patient Stations
or
Tek-CARE[®]NC200N or NC150N Nurse Call Systems
with IR151B or IR152B Patient Stations
or
Tek-CARE[®]NC110N Nurse Call System
with SF100C or SF102 Patient Stations



SF301G Call Cord

Specifications

Switch:	Momentary contact, air activated
Plug:	Stainless steel
Cord:	Vinyl tubing
Electrical:	Rated current 1/3 AMP, 24 VAC-VDC
Dimensions:	Cord: 6' long (1.83 m) Plug: 1.25" long (32 mm) Bulb: 1.5" diameter (38 mm)

Listings



UL[®] 1069 Listed and cUL[®] Listed to CSA
C22.2 N o. 205; UL[®] File Number E71599.