

SF431PSG3 Pillow Speaker IL1124 Specification Sheet

Te k-CARE400 GEN3, Tek-CARE3000 Nurse Call Systems • Rev. 4 - 2/2025

Features

- ▮ Large red nurse call button
- ▮ Three AUX call buttons for Water, Pain, and T oilet
- ▮ Four environmental control buttons
- ▮ 8-foot cord with DIN plug and security clip
- ▮ Provides pillow speaker intercom with built-in micro-
phone, environmental control, and television control
- ▮ UL[®]1069, and UL[®] 2560 Listed
- ▮ Oxygen safe
- ▮ Fully sealed for easy cleaning
- ▮ Flame-retardant (UL[®]94V-0) high-impact plastic

Description

The SF431PSG3 Pillow Speaker is used with the Tek-CARE[®]400 GEN3 IR431G3 and IR432G3 patient sta- tions, and the Tek-CARE[®]3000 IR331 and IR332 patient sta- tions.

The pillow speaker provides call placement and audio com-
munication to the Tek-CARE System and offers remote
control for hospital-grade TVs.

The pillow speaker also supplies two-way, hands-free com-
munication with the master station, and automatic silencing
of
TV audio while the intercom is in use. W ater, Pain, and Toilet
buttons, two light control buttons (typically ROOM and
READ lights), and two additional environmental controls
(ex. tem- perature up/down) are included on the pillow
speaker.

Push buttons are provided to call the nurse, plus TV
power, channel selection up/down, Closed Captioning, and
mute.

Required Components

Tek-CARE[®]400 GEN3 Nurse Call System
IR431G3 or IR432G3-series Patient Station

Tek-CARE[®]3000 Nurse Call System
IR331 or IR332 Patient Station

Optional Components

Low Voltage Controller for light and temperature control
(supplied by others)



SF431PSG3 Pillow Speaker

Listings



UL[®] 1069, UL[®] 2560, and cUL[®] listed to CSA
C22.2 No. 205; UL[®] File Number E71599.