

MODEL: R4K17V- ENHANCED SINGLE BED STATION



R4K17V

FEATURES

- One (1) DIN receptacle for patient entertainment controls and pillow speaker intercom with entertainment muting during nurse call.
- Interfaces with specialty beds via optional module
- Call-Assurance/ Monitor LED with Cancel button
- Continuous supervision of station
- Single point output for lighting control and radio stepper attachment
- Interface for full lighting control options with additional control module and pillow speakers
- Optional antimicrobial Bio-seal to cover cancel button and allow for easy disinfecting of the station
- Plug-in connectors for easy service
- Annunciates cord-out

SPECIFICATIONS

Power Required: 22 mA @15 VDC

Power Required when used with SLIM Speaker Module: 25 mA @15 VDC

Weight: 0.2 lb. (0.09 kg)

Size: W: 3.2" (8.10 cm)

H: 4.25" (10.8 cm)

D: 0.6" (1.5 cm) in back box: 0.5" (1.3 cm) off-wall

Housing and Finish: High-impact molded Cycloy 2950 UL 94-5VB; light grey faceplate

Controls/Indicators: One (1) DIN receptacle, One (1) Call-Assurance/Monitor LED, and One (1) Cancel button

Terminations: Plug-on connectors

Backbox Requirement: 1-gang, (Non-gangable) W: 2.25" (5.7 cm) H: 4" (10.2 cm) D: 3.5" (8.9 cm)

Backbox Requirement when used with SLIM Speaker Module: 3-gang, Non-gangable W: 5.5" (14 cm) H: 3.75" (9.5 cm) D: 3.5" (8.9 cm)

NOTE: If using lighting control, additional backbox space is required

Certification: UL/C-UL 1069, FCC Part 15 Class B, IEC60601-1-2, Designed for compliance to IEC60601-1 and CE Marking (MDD), FDA Registered, and RoHS Compliant

DESCRIPTION

The Responder® 4000 SLIM Enhanced Single Bed Station allows both patients and staff to place a call depending on their requirements to the assigned nurse console(s), annunciate panel(s) and or wireless devices. The Enhanced Single Bed Patient Station can interface with call cords, entertainment speakers, or specialty beds to initiate calls, control lighting, provide entertainment control, and allow intercom to the console with the optional Responder 4000 Speaker Module. When used with the NCES/NCDS Series Pillow Speakers, audio capabilities are transferred to the Pillow Speakers. When used with the Responder 4000 Feature Bed Interface module, the audio

capabilities are transferred to the bed rail with entertainment muting and the pillow speakers if used. The addition of the Responder 4000 Single Lighting Control Relay allows for full lighting control options though the buttons on the pillow speaker.

The Enhanced Single Bed Patient Station is connected to its associated corridor lamp or domeless room controller and is continually supervised for station presence. Though programming, a "Trouble" message can annunciate at the associated console(s) when the station is removed from the system.

*Architect and Engineer (A&E) Specifications available online at: customerconnection.rauland.com
Specifications subject to change without notice*

©Copyright 2009 Rauland-Borg Corporation Printed in USA Rev 03/2009