

Features

- Provides call inputs and dome light output typically for standalone devices
- Provides an audio connection in the nurse call system
- Supervised emergency and code switches
- Fully supervised wiring to the station hub
- Simplified wiring with plug-in connectors
- No visible mounting screws

Description

The IR334 Multipurpose Station provides audible, visual, and digital communications to and from the TekCARE®3000 equipment when receiving or sending signals to remotely connected devices.

The IR334 provides an output to a speaker or paging amplifier and a corridor lamp. It does not have the controls, LEDs or speakers normally associated with patient, staff or duty stations. Some practical applications for the IR334 are to provide input for code or emergency call stations in an ICU; a PA amplifier interface for the nurse call master; an alarm input point for exit/entry doors and narcotics cabinets; input for bath emergency call stations in public bath areas; or input to a remote door intercom station or access door intercom.

Required Components

NC331 Tek-CARE Hub
NC404TS and/or NC415G3 Master Stations

Optional Components

SF12x-series Peripheral Switches
IR160 Audio Station
LI484P5 Four-LED Corridor Light
PM123 Auxiliary Input Module
IH431S Surface Box



IR334 Multipurpose Station

Specifications

Dimensions:	4.80"H × 7.35"W × 2.5"D (122 mm × 187 mm × 64 mm)
Construction:	Flame-retardant (UL® 94V-0) plastic faceplate
Connections:	Plug-in 12-pin connector for up to 6 TC120 peripherals 8P8C jacks for station cabling
Housing:	2- or 3-gang electrical box and ring <i>or</i> 2- or 3-gang masonry box. Do not ground backboxes, leave floating
Wiring:	Control and power over CAT6 cable *CAT5e for retrofit

Listings



UL® 1069, UL® 2560, and cUL® listed to CSA C22.2 No. 205; UL® File Number E71599.



RoHS Compliant.