



Architects' and Engineers' Specifications

The wireless emergency code blue station shall be the part number WNC-CB and shall be ETL listed to UL1069 Edition 7 listed for "Hospital Signaling and Nurse Call Equipment." The Emergency Station shall consist of a plastic faceplate with large, blue "Code" magnetic sensor for activation and embedded transmitter. The unit shall be waterproof for use in wet areas with the addition of an included gasket. Operation of emergency signals shall be accomplished by a short pull on the nylon cord or by pulling down on the slide lever marked "Code." Once an code blue call has been placed the red LED indicator will blink until the station is reset.

System Components

WNC-SERVER-UL1069: Video Display Console
WNC-DLC: Corridor Light Power Supply
WNC-DL1W: Corridor/Zone Light
WNC-DL3WBR: Three color Corridor/Zone Light
WNC-DL4WABR: Four color Corridor/Zone Light
WNC-DS: Duty Station

WNC-CB Code Blue Station

The WNC-CB is a wireless single gang code blue pull station for use in any location where activation of a code blue signal is required. The station features an embedded transmitter and magnetic sensor and is designed to be flush mounted. The WNC-CB may also be surface mounted if desired with the addition of an enclosure and gasket. Due to the flexibility of the WNC-CB design, the station may be used as a Pull Cord Code Blue Station or a Non-Pull Cord Code Blue Station. The pull cord supplied may be used for code blue call activation, or, if desired, the slide switch may be activated. Once the call has been activated the red LED will flash until reset at the device.

Specifications

- **Dimensions:**
 Height: 4-1/2" (114mm)
 Width: 2-13/16" (72mm)
 Depth: 1-5/8" (41mm)
- Projects 11/16" (17mm) from mounting surface (including pull down lever)
- **Weight:** 5.2 oz
- **Panel:** Flame retardant GE Lexan Faceplate (UL flame rating 94V0)
- **Batteries:** Three CR123
- **Backbox:** Carlon B120A plastic single gang box
- **Flush Mount:** Cooper BB10L single gang B-Line bracket
- **Surface Mount:** Leviton 42777-1WA plastic single gang box (1.75" deep) with use of included gasket

Features

- Listed to UL1069 Ed. 7 Hospital Signaling & Nurse Call
- Simple and flexible operation
- Call-cancel at point of origin
- Flame retardant GE Lexan plastic faceplate
- Pull cord or slide switch activation
- Waterproof with the addition of gasket
- Reliable magnetic sensor
- Red LED alarm indicator