

Golden Legacy™

DHA Leg Bronzer

Australian Gold®

PREMIUM COLLECTION

FEATURES

- Leg-endary DHA Bronzer
- Flawless Finish
- Shave Minimizers
- Superior Silicones
- Caffeine
- Glycerin
- Bio-Active Bronzing Minerals
- ATO Inhibitor

PRODUCT STORY

Golden Legacy™'s Leg-endary DHA Bronzer Blend will instantly revamp your glow and gradually enhance your color for a classic golden complexion. Experience radiant, moisturized skin with the Gilded Skincare Technology. The Shave Minimizers will treat your legs to a smooth tanning base and reduce hair growth. Bio-Active Bronzing Minerals and Caffeine nourish the skin, leaving your legs unforgettably healthy and firm. The Flawless Finish will seal your tan with a soft, powder-like feel, perfect for any occasion. Let your legs leave a Golden Legacy™!

TECHNOLOGY

COLOR TECHNOLOGY

Leg-endary DHA Bronzer Blend: DHA, Caramel and Hydrolyzed Algin combine to give instant and gradual bronze results. Caramel deposits a rich hue on the top layer of the skin, while DHA continually boosts color and Hydrolyzed Algin helps increase melanin.

SKINCARE TECHNOLOGY

- **Gilded Skincare Blend:** Glycerin and a blend of Superior Silicones help provide a smooth glide while working to soothe, protect and moisturize.
- **Shave Minimizers:** Silicone, Coconut Oil and other ingredients help slow hair growth for a soft tanning base.
- **Caffeine:** Maintain everlasting bouncy skin with the antioxidant rich ingredient known to support skin tonicity.
- **Bio-Active Bronzing Minerals:** Naturally obtained from Sea Salt, these 92 essential minerals rejuvenate for visibly healthier and smoother skin.
- **Flawless Finish:** Natural Polymer leaves skin with a powder-like feel.

FRAGRANCE

Sumptuous Citrus – Strawberry, Blackberry, Apple, Peach, Jasmine, Violet, Vanilla & Amber Sandalwood.

*Refer to back of bottle for most current list of ingredients



6 FL OZ | .5 FL OZ Packette



Leg Bronzer



DHA Bronzer



ATO Inhibitor



Dermatologically Tested

HOW TO USE

This product is intended for UV tanners. Before your session, generously apply the bronzer to the legs with your hand or a mitt in a circular motion. Wash hands after application. Wait at least 2-3 hours before showering to allow for full color development. For optimal results, pair with a tan extending moisturizer to help maintain healthy, moisturized and beautiful skin. Avoid contact with eyes and delicate facial skin.

QUICK FACT

Pair Golden Legacy™ with any Premium Collection lotion for an everlasting tan!

PRICE: _____



Australian Gold.

CALIFORNIA TAN®

DESIGNER SKIN®

swedish beauty.

